
AUSWERTEN

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© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.
950M17620-011 201701001

EUROPE
Medtronic International Trading Sàrl
Route du Molliat 31
Case postale
CH-1131 Tolochenaz Switzerland
Tel: +41 (0) 21 802 70 00
Fax: +41 (0) 21 802 79 00
www.medtronic.eu

UC201902623 DEPDF

INSTALL GUIDE TO EASIER ACCESS

CareLink™ System

System set-up for administrators
and IT professionals.

ACCESS

- Users can be granted admin rights to manage user accounts. **Recommend at least 2 admins.**
- New users receive a Welcome email with their username and link to set initial password. **Password links expire after 2 hours.**
- Forgot Password functionality requires matching username and email address.

SUPPORTED ENVIRONMENTS

- Windows, Mac, or iOS
- Chrome, Internet Explorer, or Safari
- See FAQ at the bottom of <https://CareLink.Medtronic.eu> for specific versions and details

UNSUPPORTED ENVIRONMENTS

- Virtualized (e.g., Citrix)

NETWORK CONFIGURATION

- <https://CareLink.Medtronic.eu>
- IP: 144.15.107.35
- Ports 80 and 443 for regular use
- Port 7194 for one-time migration (temporary)
- **Run traceroute for 144.15.107.35 to identify any network blockages**

WORKSTATION CONFIGURATION

- Local admin rights for installation and updates
- Remove all USB devices during installation
- Close all CareLink™ applications during installation
- Disable pop-up blockers for *.CareLink.Medtronic.eu
- For silent installs, use these command line options — mode unattended — unattendedmodeui none

See next page for the steps to getting CareLink™ System online.

Medtronic

CARELINK™ SYSTEM

New web-based CareLink™ software, hosted in Medtronic-owned and controlled datacenters.

MIGRATOR


One-time copy of CareLink™ Pro patient profiles and data to the new CareLink™ System

UPLOAD

Upload current data from supported diabetes devices to the new CareLink™ System*

REPORTS

View PDF reports from supported browsers

*One time install of uploader software required

Getting new CareLink™ online is quick and easy.

Talk with your Medtronic representative today to get started.

1

Sign up

Now each employee can have his/her own account.

2

Set up

Existing data transfers to CareLink™'s secure cloud.

3

Go

Access patient data anywhere, any time.

Safety Information: CareLink™ software

CareLink™ software is intended for use as a adjunctive tool to help manage diabetes and is NOT intended to be relied upon by itself. For detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events, please consult the device manual. For further information, contact your local Medtronic representative and/or consult the User Guides at medtronicdiabetes.com/support/download-library/user-guides.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

EUROPE

Medtronic International Trading Sàrl
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Switzerland
Tel: +41 (0) 21 802 70 00
Fax: +41 (0) 21 802 79 00

www.medtronic.eu

FREQUENTLY ASKED QUESTIONS CARELINK[®] SYSTEM

Medtronic

Who can use this system?

- Healthcare professionals with a registered clinic and valid user account can log in to this system.

Why should I use this system?

- The CareLink software helps healthcare professionals make more informed diabetes therapy decisions. CareLink reports combine all blood glucose, sensor glucose, and insulin pump data into a complete picture, with intelligent analysis and pattern recognition that provide actionable insights. Furthermore, the new CareLink software provides easier access from more devices (e.g. Mac and iPad) and an updated design for a better user experience.

How can I register to use this system?

- Contact your local Medtronic sales representative for your unique, single-use registration code.
- During registration, the first registered user will be assigned clinic admin privileges
 - Clinic admin privileges can create new user accounts, manage existing accounts, and manage clinic default settings.
 - Recommend at least two clinic admins at all times for contingencies.

My IT and/or Legal departments have questions. Who can they contact?

- We have created resources and processes to assist with IT and Legal questions. Please contact our Medtronic 24-hour helpline for more information. The Medtronic 24-hour helpline is available 24/7 at (800)646-4633.

What are the computer system requirements to use this system?

- Operating Systems
 - Microsoft[®] Windows[®] 7, 8.1, 10
 - Mac OS[®] 10.10 Yosemite, 10.11 El Capitan, 10.12 Sierra
 - iOS 9.3
- Browsers
 - Internet Explorer[®] 11
 - Safari[®] 9, 10
 - Chrome[™]

What devices can I upload?

- Meters
 - Ascensia
 - CONTOUR®NEXT LINK
 - CONTOUR®NEXT LINK 2.4
 - BREEZE®
 - CONTOUR®NEXT
 - CONTOUR®NEXT EZ
 - CONTOUR®NEXT USB
 - CONTOUR®
 - CONTOUR®LINK
 - CONTOUR®USB
 - LifeScan
 - OneTouch® Ultra
 - OneTouch® Ultra2
 - OneTouch® UltraLink
 - OneTouch® UltraSmart
 - OneTouch® UltraMini/UltraEasy
 - Abbott
 - Freestyle®
 - Freestyle® Lite
 - Roche
 - Accu-chek® Active
 - Accu-chek® Compact
- Insulin Pumps
 - MiniMed® Paradigm® REAL-Time Revel™ system
 - MiniMed® 530G system
 - MiniMed® 630G system
 - MiniMed® 670G system
 - MiniMed® Paradigm® 515/715, 522/722, 523/723, 551/751 systems
- Continuous Glucose Monitoring
 - Guardian® REAL-Time system

What do I need to upload a device?

- Communications devices
 - CareLink USB device
 - NOTE: The white CareLink USB device used for Paradigm pump systems (e.g. MiniMed 530G) is not supported on Mac. To upload these pumps on a Mac, please use the Ascensia CONTOUR®NEXT LINK meter as the communications device.
 - Ascensia CONTOUR®NEXT LINK for MiniMed 530G system

- Ascensia CONTOUR®NEXT LINK 2.4 for MiniMed 630G system or MiniMed 670G system
- Manufacturer supplied meter cables
- Software
 - CareLink software needs the uploader application installed in order to add devices and upload patient data. The first time you click the upload button, you will be presented with information to download and install the uploader application. Please follow the on-screen instructions to complete the installation.
 - The Java applet is not required or used by this CareLink software to upload.
- Troubleshooting
 - To install the uploader application as an administrator on a Mac, you may need to hold the "Ctrl" button and click the installation file. Then, you can select the "Open" option to continue with the installation.
 - Nearing the end of the installation, you will be prompted with a message to disconnect any medical or USB devices. Please ensure that you do not have any meters or CareLink USB devices plugged in and click "OK" to continue the installation.
 - When you click the upload button, your computer may send a message asking if you would like to launch this application. Click to "Allow" this action to avoid further pop-up blocking.
- Linking
 - If your patients upload to CareLink Personal software, you can automatically download the latest data to your CareLink software so that your patients can upload at their convenience and you do not need to spend time uploading in your office. You can set up a link under the patient profile section of your CareLink software. Your patient's CareLink Personal username and password must be entered when setting up this link.

After I upload, why doesn't my data calendar show any new or updated data?

- After completing the upload, to refresh the data calendar, please click on the "Refresh" button on the right side of the screen.
- If you need additional assistance on this screen, click on the "i" button found to the right of the Reporting Period label.

What happens to all my existing historical CareLink Pro data and patients?

- Data migration helps you get started by copying your patient profiles and data from your previous CareLink Pro desktop software to this new CareLink online software. A guide on *How to Migrate* is accessible via a link at the top of the Clinic Dashboard screen. For more information on migration, please refer to the Help section within the application.

Where will my patients' data be stored?

- The new CareLink online software uses the latest technology to store and safeguard your patient's data on a server that Medtronic maintains.

What do I do if I forget my username?

- Contact your clinic administrator.

What do I do if I forget my password?

- A "Forgot your password?" link is available on the home page. After successfully submitting your request, an email will be sent to the email address on file with a link to reset your password.
- If you have issues with the "Forgot your password?" process, you may contact your clinic administrator to send a password reset email.

I tried to "Generate" reports, but nothing happened?

- You may need to disable your browser's pop-up blocker. Please refer to each specific browser's pop-up blocker instructions for proper operation. Some suggested resources are below.
- For Internet Explorer: <https://support.microsoft.com/en-us/help/17479/windows-internet-explorer-11-change-security-privacy-settings>
- For Safari: https://support.apple.com/kb/PH21485?locale=en_US
- For Chrome: <https://support.google.com/chrome/answer/95472?co=GENIE.Platform%3DDesktop&hl=en>
- If you need additional assistance on this screen, click on the "i" button found to the right of the Reporting Period label.

Why do I need a clinic administrator?

- Clinic administrators can assign admin privileges, add new user accounts, send password reset emails, and perform other admin-only functions.

How many clinic administrators do you need?

- Recommend at least two clinic admins at all times for contingencies.

CONTOUR and BREEZE are registered trademarks of Ascensia Diabetes Care Inc. OneTouch is a registered trademark of LifeScan, Inc, a Johnson & Johnson company. Freestyle is a registered trademark of Abbott Diabetes Care, Inc. Accu-chek is a registered trademark of Roche Diabetes Care, Inc. Microsoft and Internet Explorer are registered trademarks of Microsoft, Inc. Mac OS and Safari are registered trademarks of Apple, Inc. Chrome is a trademark of Google, Inc. CareLink, MiniMed, Paradigm, and Guardian are registered trademarks, and Revel is a trademark of Medtronic MiniMed, Inc. 940M16211-011 20170209 ©2017 Medtronic MiniMed, Inc. All rights reserved.